

# 3M 8541 Brown pad

## For use with Doodlebug™ pad holder



Technical data sheet

TDS 135

Date: May 2024

### Description:

The product 3M 8541 Brown Doodlebug™ pad is a thick brown stripping pad, made of high-quality synthetic fibre forming an open textured non-woven web, having abrasive particles evenly dispersed throughout and bonded to the web with a durable resin.

### Product Advantages:

- Excellent stripping effect
- Consistent performance
- Flexible and able to conform to uneven surface
- Long lasting.
- Easy to rinse/wash after use.
- Suitable for dry or wet stripping operations
- Excellent resistance to stripping solutions and 3M Trouble Shooter

### Where and When to Use:

Use the pad attached to a Doodlebug™ pad holder to remove marks, scuffs and stains from non-fragile floors, corners, baseboards, walls and other surfaces.

Maximum use/washing temperature: 40°C.

### Labelling:

Each case shall be labelled to identify the content by product name, size and manufacturer.

### Packaging:

30 pads per case (5 pads in multilingual printed pack, 6 intermediates per pack).

### Maintenance:

Remove from holder and rinse the pad after each use under tap water.

Brush the pad with a nylon brush to remove deep embedded soil if necessary.

### Disposal information:

Discard into appropriate container for incineration, if available.

### Product Characteristics:

<i>Colour:</i>	Brown
<i>Thickness:</i>	28 mm (*)
<i>Size:</i>	254 mm x 117 mm
<i>Weight:</i>	2300 g/m <sup>2</sup> (*)
<i>Composition:</i>	Synthetic fibres, minerals, polymer resin.

### Aggressiveness:

2,4 g removed (\*)

(Schiefer test on acrylic plate material, after 5000 revolutions with a new sample pad)

(\*) Typical average values

### Origin:

Made in Thailand in an ISO 9001:2015 and ISO14001:2015 3M certified plant.