













Mini Alexander

Bracket replacement code

U = UPPER

Tooth	Prescription	Torque	Angulation	Height	Width	Base	ID Color	Slot .018x.028		Slot .022x.028	
								Right	Left	Right	Left
Central	Mini Alexander	+15°	5°	1.6	3.5	-		210-001	210-002	210-101	210-102
Lateral	Mini Alexander	+8°	8°	1.7	2.8	-		210-011	210-012	210-111	210-112
Cuspids with hooks	Mini Alexander	-2°	11°	1.6	6.0	-		210-021	210-022	210-121	210-122
Bicuspid with hooks	Mini Alexander	-7°	0°	1.6	5.2	-		210-031	210-031	210-131	210-131
Bicuspid without hooks	Mini Alexander	-7°	0°	1.6	5.2	-		210-033	210-033	210-133	210-133

L = LOWER

Tooth	Prescription	Torque	Angulation	Height	Width	Base	ID Color	Slot .018x.028		Slot .022x.028	
								Right	Left	Right	Left
Central Incisors without hooks	Mini Alexander	-4°	0°	2.0	4.4	-		210-045	210-045	210-145	210-145
Lateral Incisors	Mini Alexander	-4°	0°	2.0	4.4	-		210-046	210-046	210-146	210-146
Cuspids with hooks	Mini Alexander	-7°	6°	1.5	6.0	-		210-051	210-052	210-151	210-152
1 st Bicuspid with hooks	Mini Alexander	-11°	0°	1.7	5.2	-		210-061	210-061	210-161	210-161
1 st Bicuspid without hooks	Mini Alexander	-11°	0°	1.7	5.2	-		210-063	210-063	210-163	210-163
2 nd Bicuspid with hooks	Mini Alexander	-17°	0°	1.6	5.2	-		210-071	210-071	210-171	210-171
2 nd Bicuspid without hooks	Mini Alexander	-17°	0°	1.6	5.2	-		210-073	210-073	210-173	210-173