



3M™ Red Dot™ Multi-Purpose Monitoring Electrodes with Sticky Gel 2560 Electrode

NPC Code	FDK001	
Indication for Use	Long Term Monitoring	
Dimensions	Electrode Size: Skin Contact Size: Adhesive Area: Height excluding connector:	4.06cm x 3.45cm 4.06cm x 3.45cm 10.51 sq cm 0.19cm
Electrode Materials	Backing Material: Backing Material Adhesive: Connector: Release Liner: Sensor Material:	Foam 3M™ Durapore™ Adhesive Stainless Steel Snap Si coated paper Silver/silver-chloride coated plastic
Sensor	Gel system: Gel area: Sensor area:	Sticky gel 3.48 sq cm 3.48 sq cm
Lifetime	Recommended max application time: Sealed pouch:	5 days 2 years
X-Ray and MRI	X-Ray: MRI:	NO NO
Environmental Issues	PVC-free electrode: Latex-free electrode: PVC-free packaging	YES YES YES
Packaging Quantities	Pouch:	50

3M Health Care Limited

3M House
Morley Street
Loughborough
Leicestershire LE11 1EP
Telephone (01509) 611611
Facsimile (01509) 237288
www.3mhealthcare.co.uk

3M Ireland

3M House, Adelphi Centre
Dun Laoghaire
Co Dublin
Ireland
Telephone (01) 280 3555
Facsimile (01) 280 3509

3M and Red Dot are trademarks
of the 3M Company.
© 3M Health Care 2008.

