

Fasteners, Closures & Components

Product Selection Guide

Low Profile Fasteners/Closures

Properties	Plain Back, No Adhesive		With Adhesive				
	7330	7334	7331	7335	7337	7336	
1. Closure Type	Loop	Hook	Loop	Hook	Loop	Hook	
2. Sterilization	G, E	E	G, E	E	G, E	E	
3. Color	White	White	White	White	White	White	
4. Heat Sealable	Yes	Yes	N/A	N/A	N/A	N/A	
5. Dynamic Shear Closure 0.5 inch ² /3.25mm ²	290 oz (8.3kg)	290 oz (8.3kg)	290 oz (8.3kg)	290 oz (8.3kg)	290 oz (8.3kg)	290 oz (8.3kg)	
6. 135 deg Peel Closure	12.5 oz/inch wide (360g/25mm wide)		12.5 oz/inch wide (360g/25mm wide)		12.5 oz/inch wide (360g/25mm wide)		
7. Closure Caliper	mil	9.8	18.5	14	22.7	10	19.0
	millimeters	0.25	0.5	0.4	0.6	0.25	0.5
8. Adhesive Description	N/A	N/A	High Tack, Syn. Rubber-Based	High Tack, Syn. Rubber-Based	High Tack, Syn. Rubber-Based	High Tack, Syn. Rubber-Based	
9. Adhesion to	oz/inch	N/A	N/A	36 (LDPE)	36 (LDPE)	90 (S.Steel)	80 (S.Steel)
	kg/25mm	N/A	N/A	1.0 (LDPE)	1.0 (LDPE)	2.6 (S.Steel)	2.3 (S.Steel)
10. Release Liner Caliper	mil	N/A	N/A	3.2	3.2	3.5	3.5
	millimeters	N/A	N/A	0.08	0.08	0.09	0.09
11. Width	inches	0.5 – 30	0.5 – 28	0.5 – 28	0.5 – 26	0.5 – 28	0.5 – 26
	centimeters	1.3 – 76	1.3 – 71	1.3 – 71	1.3 – 66	1.3 – 71	1.3 – 66
12. Length	yards	< 1.5"	175	150	150	150	150
		> 1.5"	450	250	145	145	145
	meters	< 3.8cm	160	137	137	137	137
		> 3.8cm	411.5	228.6	132	132	132
13. Approx. Roll Diameter	inches	< 1.5"	10	10	10	10	10
		> 1.5"	15	15	15	15	15
	centimeters	< 3.8cm	25	25	25	25	25
		> 3.8cm	38	38	38	38	38

Products 7334 and 7335 may not be gamma-stable at all customers' levels of exposure (see page 69A). Customers will need to evaluate their device construction requirements using the replacement. Please contact the 3M Center Medical Specialties group for more information and samples.

Table Revised November 2012

Tube & Cord Organizer

Product Number	Color	Width in. (cm)	Total Length in. (cm)	Organizers/ Carton	Typical Adhesion to Stainless Steel oz (kg)
1115NS	White	1.5 (3.8)	5.25 (13.3)	500	60 (1.7)

This is a tape-based recloseable organizer for use on surgical gowns and drapes or equipment. The product is supplied in bulk packaging, can be EtO sterilized for use in the operating room and included in procedure kits and packs.