

3M Food Safety
3M™ Prepared Enrichment Media

Collect ▶ Transport ▶ Prepare ▶ Analyze ▶

Experience confidence

The complete 3M sample handling system

Prepare.

3M Food Safety offers a full line of prepared enrichments and diluents that provide a cost-effective alternative to making and sterilizing your own media. Standardize your testing results with our wide variety of prepared media, available in convenient flip-top vials, pouches or bottles.

- Convenient, ready-to-use media for consistent, reliable results
- Cost-effective alternative to making and sterilizing media
- Lowers operational labor costs
- Reduces errors in mixing and preparing dilutions
- Wide range of broths and buffers

The 3M logo is located in the bottom right corner of the page. It consists of the letters '3M' in a bold, red, sans-serif font.

Prepared **enrichments** and **diluents** for **standardized test results** — **test after test.**



Collect ▶

Transport ▶

Prepare ▲

Analyze ▶

3M™ Flip-Top Dilution Bottles



- Flip-top lids are easy to open
- Hinged lid stays open for easy sample addition
- 45mm bottle opening makes it easy to pipet into and out of the bottle
- Choose 90mL or 99mL fill volumes
- Available with Butterfield's Buffer, Peptone Water, Buffered Peptone Water or Phosphate Buffer for increased flexibility
- Two-year shelf life from date of manufacture

3M™ Mini Flip-Top Vials



- Pre-filled vials contain 9mL or 10mL of diluents
- Miniaturized dilution bottles are less expensive to ship and conserve valuable lab space
- Flip-top lids are easy to open
- Hinged lid stays open for easy sample addition
- Available with a variety of broths and buffers

3M™ Dilution and Enrichment Pouches



- Just add sample, homogenize to mix, then place directly into the incubator
- Pre-measured with 99mL, 225mL or 450mL of broths or buffers
- Strong, multi-laminate pouch resists punctures
- Sturdy pouch stands alone
- Wire tabs keep bag open for easy sample addition

Ordering Information

	Catalog No.	Description	Qty.
Flip-Top	FTBFD90	90mL Butterfield's Buffer	66
	FTBFD99	99mL Butterfield's Buffer	66
	FTPW9066	90mL Peptone Water	66
	FTPW99	99mL Peptone Water	66
	FTPHB9066	90mL Phosphate Buffer	66
	FTPHB9966	99mL Phosphate Buffer	66
	FTBW9066	90mL Buffered Peptone Water	66
	FTBPWIS090	90mL Buffered Peptone Water (ISO)	66
Mini Flip-Top	BPPFV9BPW	9mL Buffered Peptone Water Broth	100
	BPPFV9BFD	9mL Butterfield's Buffer	100
	BPPFV9LB	9mL Lethen Broth	100
	BPPFV9MRD	9mL Maximum Recovery Diluent	100
	BPPFV9MRS	9mL MRS Broth	100
	BPPFV10NB	10mL Neutralizing Buffer	100
	BPPFV10BPW	10mL Buffered Peptone Water Broth	100
	BPPFV10DE	10mL D/E Broth	100
BPPFV10LB	10mL Lethen Broth	100	
Dilution	QDBFD99	99mL Butterfield's Buffer	50
	QDBFD225	225mL Butterfield's Buffer	32
	QDBFD225F	225mL Butterfield's Buffer w/Filter	32
	QDBFD450	450mL Butterfield's Buffer	16
Enrichment	QEBPW225	225mL Buffered Peptone Water Broth	32
	QEBPW225F	225mL Buffered Peptone Water Broth w/Filter	32
	QEDFB225	225mL Demi Fraser Broth	32
	QELAC225	225mL Lactose Broth	32
	QELAC225F	225mL Lactose Broth w/Filter	32
	QEMBPW225	225mL Modified Buffered Peptone Water Broth	32
	QETSB225	225mL Tryptic Soy Broth	32
	QEUV225	225mL UVM	32
QEBPWIS225	225mL Buffered Peptone Water (ISO)	32	
QEBPWIS225F	225mL Buffered Peptone Water (ISO) w/Filter	32	



3M Food Safety
3M Center
Building 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553
www.3M.com/foodsafety

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-364-3577

For more information, contact your local 3M Food Safety representative or visit 3M.com/foodsafety.

Please recycle. Printed in U.S.A.
© 3M 2013. All rights reserved.
70-2011-5013-6

3M is a trademark of 3M.
Used under license in Canada.