

## Adhesives

3M™ Scotchbond™ Universal Adhesive

**3+1**

Buy 3 of the Same, Get 1 of the Same Free!



## CAD/CAM

3M™ Chairside Zirconia

**4+1**

Buy Any 4, Get 1 Free!\*



## Cements

3M™ RelyX™ Cements

3M™ RelyX™ Luting Plus Resin Modified Cement

3M™ RelyX™ Ultimate Adhesive Resin Cement

3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement

**4+1**



Buy Any 4, Get 1 Free!\*

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement

**4+1**

Aplicap™ / Maxicap™

Buy Any 4, Get 1 Free!\*



## Composites

3M™ Filtek™ Restoratives

3M™ Filtek™ Bulk Fill Flowable Restorative

3M™ Filtek™ One Bulk Fill Restorative

3M™ Filtek™ Supreme Ultra Flowable Restorative

3M™ Filtek™ Supreme Ultra Universal Restorative

3M™ Filtek™ Universal Restorative

**3+1**

Buy Any 3 Refills, Get 1 Free!\*



## Impressioning

3M™ Paradigm™ VPS Impression Material

**4+1**

Buy Any 4 Refills, Get 1 Free!



## Impressioning

3M™ Penta™ Impression Material

3M™ Impregum™ Penta™ Super Quick Polyether Impression Material

3M™ Impregum™ Penta™ Soft Polyether Impression Material

3M™ Impregum™ Penta™ Polyether Impression Material

3M™ Imprint™ 4 Penta™ VPS Impression Material

3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material

**4+1** Penta™

Buy Any 4, Get 1 Free!\*



3M™ Garant™ Impression Material

3M™ Impregum™ Super Quick Polyether Impression Material

3M™ Impregum™ Soft Polyether Impression Material

3M™ Imprint™ 4 VPS Impression Material

3M™ Imprint™ 4 Preliminary VPS Impression Material

**4+1** Garant™

Buy Any 4, Get 1 Free!\*



3M™ Imprint™ 4 Bite VPS Registration Material

**4+1**

Buy 4 Refills, Get 1 Free!



## Preventive

3M™ Clinpro™ 5000 Anti-Cavity Toothpaste

**4+1**

Buy Any 4 Cases, Get 1 Case Free!



3M™ Vanish™ White Varnish

3M™ Fast Release Varnish

**4+1**

Buy Any 4 Same Count Packs, Get 1 of the Same Free!



3M™ Xerostomia Relief Spray

**15+2**

Buy 15 Boxes, Get 2 Free!



**3M Oral Care Special Offers Bonus Product Redemption.**

You can redeem for your bonus goods by email. Refer to the instructions below. Proof of purchase on one 3M Oral Care Authorized Distributor invoice with order dated January 1 – March 31, 2020 is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). Redeem by April 30, 2020. Offers cannot be combined with other special offers (e.g., DentaCheques). Offer void if purchased product is returned. Bonus goods can only be returned within 30 days of ship date.

**Redemption by email: 3MOralCarePromos@mmm.com**

Complete this coupon and email it with your invoice by April 30, 2020. Please provide complete shipping information, including business phone.

BUSINESS PHONE		FAX
DOCTOR NAME		
BUSINESS ADDRESS		
CITY	STATE	ZIP
EMAIL ADDRESS		

**A. 3M™ Scotchbond™ Universal Adhesive**  
*Buy 3 of the Same, Get 1 of the Same Free*  
 Vial Intro Kit (41254)  
 Contains: 1 – 5ml adhesive vial, 50 mixing wells, 50 microbrushes, 3ml etchant syringe, 25 dispensing tips, step-by-step cards  
**Vial Refill (41258)**  
 Contains: 1 – 5ml adhesive vial  
**Unit Dose Packs**  
 Contains: unit dose applicators, 3ml etchant syringe, 25 dispensing tips, step-by-step cards  
**40-Pack (41255), 100-Pack (41256), 200-Pack (41257)**

**B. 3M™ Chairside Zirconia** *Buy Any 4, Get 1 Free\**  
**Intro Kit (69956)**  
 Contains: 1 crown block A1, 2 crown blocks A2, 1 crown block A3, 1 bridge block A2  
**Crown Block Refills**  
 Contains: 5 – 20mm blocks  
**Bleach (69940), A1 (69941), A2 (69942), A3 (69943), A3.5 (69944), B1 (69945), C1 (69946), D2 (69947)**  
**Bridge Block Refills**  
 Contains: 3 – 39mm blocks  
**Bleach (69948), A1 (69949), A2 (69950), A3 (69951), A3.5 (69952), B1 (69953), C1 (69954), D2 (69955)**

**C. 3M™ RelyX™ Cements** *Buy Any 4, Get 1 Free\**  
**3M™ RelyX™ Luting Plus Resin Modified Cement**  
**Automix Value Pack (3535)**  
 Contains: 3 – 8.5g syringes, 36 mixing tips  
**Automix Trial Kit (3535TK)**  
 Contains: 1 – 8.5g syringe, 12 mixing tips  
**Clicker™ Refill (3525)**  
 Contains: 2 – 11g dispensers  
**Clicker™ Trial Kit (3525TK)**  
 Contains: 1 – 11g dispenser  
**3M™ RelyX™ Ultimate Adhesive Resin Cement**  
**Refills**  
 Contains: 1 – 8.5g syringe cement, 10 mixing tips (regular), 5 mixing tips (wide), 5 intra-oral tips  
**A1 (56887), A3O Opaque (56888), B0.5 Bleach (56889), Translucent (56890)**  
**Intro Kit (56891)**  
 Contains: 1 – 5ml vial Scotchbond™ Universal Adhesive, 1 – 3ml etchant syringe, 25 dispensing tips for etchant syringe, 50 mixing wells, 50 microbrushes, 1 – 8.5g cement syringe each of shades A3O, B0.5, A1 and Translucent, 30 mixing tips (regular), 10 mixing tips (wide) + 10 endo tips, 10 mixing tips (wide) + 10 intra-oral tips  
**Trial Kits**  
 Contains: 1 – 5ml vial Scotchbond™ Universal Adhesive, 1 – 3ml etchant syringe, 25 dispensing tips for etchant syringe, 50 microbrushes, 1 – 8.5g cement syringe, 10 mixing tips (regular), 5 mixing tips (wide) + 5 intra-oral tips  
**A1 (56894), Translucent (56892)**

**3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement**  
**Automix Refills**  
 Contains: 1 – 8.5g syringe, 10 mixing tips (regular), 5 mixing tips (wide), 5 endo tips  
**A2 Universal (56846), A3 Opaque (56847), Translucent (56848)**  
**Automix Value Packs**  
 Contains: 3 – 8.5g syringes, 30 mixing tips (regular), 20 mixing tips (wide), 15 endo tips, 5 intra-oral tips  
**A2 Universal (56849), A3 Opaque (56857), Translucent (56858)**  
**Automix Starter Kit (56850)**  
 Contains: 1 – 8.5g syringe A2 universal, 7 mixing tips (regular), 15 mixing tips (wide), 15 endo tips, 5 intra-oral tips, RelyX™ Fiber Posts 5 of each size 0, 1, 2 and 3, drills sizes 0, 1, 2, 3 and Universal  
**Clicker™ Refills**  
 Contains: 1 – 11g dispenser, mixing pad  
**A2 Universal (56875), A3 Opaque (56876), Translucent (56874)**

**D. 3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement** *Buy Any 4, Get 1 Free\**  
**Aplicap™ Refills – 50 capsules**  
**Assorted (56815), Translucent (56816), A3 Opaque (56817), A2 Universal (56818)**  
**Aplicap™ Refills – 20 capsules**  
**A1 (56819), White Opaque (56820)**  
**Aplicap™ Intro Kit (56814)**  
 Contains: 50 assorted capsules, activator, applicator  
**Aplicap™ Trial Kit (56829)**  
 Contains: 20 assorted capsules, activator, applicator  
**Maxicap™ Refills – 20 capsules**  
**Assorted (56831), Translucent (56832), A3 Opaque (56833), A2 Universal (56834)**  
**Maxicap™ Intro Kit (56830)**  
 Contains: 20 assorted capsules, activator, applicator

**E. 3M™ Filtek™ Restorative Refills**  
*Buy Any 3 Refills, Get 1 Free\**  
 Specify refills:  
**3M™ Filtek™ Bulk Fill Flowable Restorative**  
 Syringe 2 – 2g (4862 series)  
 Capsule 15 – 0.2g (4861 series)  
**3M™ Filtek™ One Bulk Fill Restorative**  
 Syringe 1 – 4g (4866 series)  
 Capsule 20 – 0.2g (4867 series)  
**3M™ Filtek™ Supreme Ultra Flowable Restorative**  
 Syringe 2 – 2g (6032 series)  
 Capsule 20 – 0.2g (6033 series)  
**3M™ Filtek™ Supreme Ultra Universal Restorative**  
 Syringe 1 – 4g (6028 series)  
 Capsule 20 – 0.2g (6029 series)  
 Capsule 10 – 0.2g (6029 series)  
**3M™ Filtek™ Universal Restorative**  
 Syringe 1 – 4g (6555 series)  
 Capsule 20 – 0.2g (6550 series)  
 Capsule 10 – 0.2g (6550 series)

**F. 3M™ Paradigm™ VPS Impression Material**  
*Buy Any 4 Refills, Get 1 Free*  
 Contains: 2 – 50ml cartridges  
**Heavy Body (5313), Heavy Body Fast Set (5513), Light Body (5311), Light Body Fast Set (5511), Monophase (5312), Regular Body (5315), Regular Body Fast (5515)**

**3M Oral Care Products**  
**Find a U.S. Authorized**  
**3M Oral Care Distributor at:**  
**3M.com/BuyDental**

**G. 3M™ Penta™ Impression Material**  
*Buy Any 4, Get 1 Free\**  
**3M™ Impregum™ Penta™ Super Quick Polyether Impression Material**  
**Refills**  
 HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst Super Quick Heavy Body (69406), Super Quick Medium Body (69385)  
**3M™ Impregum™ Penta™ Soft Polyether Impression Material**  
**Refills**  
 HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst LB Contains: 1 – 300ml base, 1 – 60ml catalyst  
**Quick Step Medium Body (31765), Quick Step Heavy Body (31769), Light Body (31746), Medium Body (31731), Heavy Body (31741)**  
**Intro Kits**  
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta™ Cartridge, 1 Penta™ Elastomer Syringe, 10 Penta™ Mixing Tips, 1 polyether tray adhesive  
 HB/LB Also Contains: 2 – 50ml cartridges LB, 5 polyether mixing tips, 5 polyether intraoral tips, 5 intra-oral syringes  
**Heavy Body/Light Body (P31739) + 1 cartridge dispenser, Quick Step Heavy Body/Light Body (P31779), Medium Body (P31735), Quick Step Medium Body (P31771)**  
**3M™ Impregum™ Penta™ Polyether Impression Material**  
**Refill Medium Body (31642)**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Intro Kit Medium Body (P31685)**  
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta™ Cartridge, 1 Penta™ Elastomer Syringe, 10 Penta™ Mixing Tips, 1 polyether tray adhesive  
**3M™ Imprint™ 4 Penta™ VPS Impression Material**  
**Refills**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Putty (71486), Heavy Body (71484), Super Quick Heavy Body (71485)**  
**Intro Kits**  
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 1 Penta™ Cartridge, 10 Penta™ Mixing Tips, 5 intra-oral syringes, 5 mixing tips, 5 intra-oral tips, 1 VPS tray adhesive  
**Putty (71536), Heavy Body (71494), Super Quick Heavy Body (71495)**  
**3M™ Imprint™ 4 Preliminary Penta™ VPS Preliminary Impression Material**  
**Refills**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Regular (71521), Super Quick (71522)**  
**Value Packs**  
 Contains: 6 – 300ml base, 6 – 60ml catalyst  
**Regular (71523), Super Quick (71524)**

**H. 3M™ Garant™ Impression Material** *Buy Any 4, Get 1 Free\**  
**3M™ Impregum™ Super Quick Polyether Impression Material**  
**Refill Super Quick Light Body (69379)**  
 LB Contains: 4 – 50ml cartridges, 5 intra-oral syringes, 5 mixing tips  
**3M™ Impregum™ Soft Polyether Impression Material**  
**Refills**  
 Tray Contains: 4 – 50ml cartridges, 10 mixing tips  
 LB Contains: 4 – 50ml cartridges, 10 mixing tips, 5 intra-oral syringes  
**Tray (31785), Quick Step Tray (31788), Light Body (31752), Quick Step Light Body (31767)**  
**Intro Kits**  
 Contains: 3 – 50ml tray cartridges, 1 – 50ml LB cartridge, 10 polyether tray mixing tips, 5 polyether mixing tips, 5 polyether intraoral tips, 5 intra-oral syringes, 1 polyether tray adhesive  
**Tray (31786), Quick Step Tray (31789)**  
**3M™ Imprint™ 4 VPS Impression Material**  
**Refills**  
 HB Contains: 4 – 50ml cartridges, 10 mixing tips  
 RB/LB Contains: 4 – 50ml cartridges, 5 mixing tips, 5 intra-oral syringes  
**Heavy Body (71492), Light Body (71488), Regular Body (71487), Super Quick Heavy Body (71493), Super Quick Light Body (71490), Super Quick Regular Body (71489)**  
**Intro Kits**  
 Contains: 2 – 50ml HB cartridges, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 10 mixing tips, 10 intra-oral syringes, 1 VPS tray adhesive  
**Heavy Body (71534), Super Quick Heavy Body (71535)**  
**Economy Packs**  
 Contains: 10 – 50ml cartridges  
**Heavy Body (71531), Super Quick Heavy Body (71532)**  
**3M™ Imprint™ 4 Preliminary VPS Preliminary Impression Material**  
**Economy Packs**  
 Contains: 10 – 50ml cartridges  
**Regular (71539), Super Quick (71540)**  
**Bulk Packs**  
 Contains: 30 – 50ml cartridges  
**Regular (71548), Super Quick (71549)**

**I. 3M™ Imprint™ 4 Bite VPS Bite Registration Material Refill (71529)** *Buy 4 Refills, Get 1 Free*  
 Contains: 2 – 50ml cartridges, 5 mixing tips, 5 intra-oral tips bite

**J. 3M™ Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste**  
*Buy Any 4 Cases, Get 1 Case Free*  
 Contains: 24 – 4 oz. tubes per case  
**Vanilla Mint (12115), Spearmint (12115SM), Bubble Gum (12115BG)**

**K. Buy Any 4 Same Count Varnish Packs, Get 1 of the Same Free**  
**3M™ Vanish™ 5% Sodium Fluoride White Varnish**  
 50-Pack: Cherry (12149C), Mint (12149M), Melon (12149L)  
 100-Pack: Cherry (12150C), Mint (12150M), Melon (12150L)  
 1,000-Pack: Cherry (12151C), Mint (12151M), Melon (12151L), Cherry/Melon Mix (12151CL), Cherry/Melon/Mint Mix (12151X)  
**3M™ Fast Release Varnish 5% Sodium Fluoride**  
 50-Pack: Mint (12157M), Strawberry (12157S), Caramel (12157C)  
 200-Pack: Mint (12158M), Strawberry (12158S), Caramel (12158C)

**L. 3M™ Xerostomia Relief Spray (12137)**  
*Buy 15 Boxes, Get 2 Free*  
 Contains: 2 – 10ml spray canisters

**A. Get one of the Same Scotchbond Universal Adhesive kits, refills or packs purchased**  
 Vial: \_\_\_Intro Kit (41254) \_\_\_Refill (41258)  
 Unit Dose: \_\_\_40-Pack (41255) \_\_\_100-Pack (41256) \_\_\_200-Pack (41257)

**B. Get one Chairside Zirconia\***  
 Intro Kit: \_\_\_(69956)  
 Crown Block Refill: \_\_\_Bleach (69940) \_\_\_A1 (69941) \_\_\_A2 (69942) \_\_\_A3 (69943) \_\_\_A3.5 (69944) \_\_\_B1 (69945) \_\_\_C1 (69946) \_\_\_D2 (69947)  
 Bridge Block Refill: \_\_\_Bleach (69948) \_\_\_A1 (69949) \_\_\_A2 (69950) \_\_\_A3 (69951) \_\_\_A3.5 (69952) \_\_\_B1 (69953) \_\_\_C1 (69954) \_\_\_D2 (69955)

**C. Get one RelyX Cement\***  
**RelyX Luting Plus Resin Modified Cement**  
 Automix: \_\_\_Value Pack (3535) \_\_\_Trial Kit (3535TK)  
 Clicker: \_\_\_Refill (3525) \_\_\_Trial Kit (3525TK)  
**RelyX Ultimate Adhesive Cement**  
 Refill: \_\_\_A1 (56897) \_\_\_A3O Opaque (56888) \_\_\_B0.5 Bleach (56889) \_\_\_Translucent (56890)  
**Intro Kit: (56891)**  
 Trial Kit: \_\_\_A1 (56894) \_\_\_Translucent (56892)  
**RelyX Unicem 2 Cement**  
 Automix Refill: \_\_\_A2 Universal (56846) \_\_\_A3 Opaque (56847) \_\_\_Translucent (56848)\*\*  
 Automix Value Pack: \_\_\_A2 Universal (56849) \_\_\_A3 Opaque (56857) \_\_\_Translucent (56858)  
 Automix Starter Kit: \_\_\_(56850)  
 Clicker Refill: \_\_\_A2 Universal (56875) \_\_\_A3 Opaque (56876) \_\_\_Translucent (56874)

**D. Get one RelyX Unicem Self-Adhesive Universal Resin Cement\***  
**Aplicap**  
 Refill: \_\_\_Assorted (56815) \_\_\_Translucent (56816)\*\* \_\_\_A3 Opaque (56817) \_\_\_A2 Universal (56818) \_\_\_A1 (56819) \_\_\_White Opaque (56820)  
**Intro Kit: (56814) Trial Kit: (56829)**  
**Maxicap**  
 Refill: \_\_\_Assorted (56831) \_\_\_Translucent (56832) \_\_\_A3 Opaque (56833) \_\_\_A2 Universal (56834)  
**Intro Kit: (56830)**

**E. Get one Filtek Restorative refill\***  
**Filtek Bulk Fill Flowable Restorative**  
 Choose \_\_\_syringe (4862) or \_\_\_capsule (4861):  
 \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_U  
**Filtek One Bulk Fill Restorative**  
 Choose \_\_\_syringe (4866) or \_\_\_capsule (4867):  
 \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_B1 \_\_\_C2  
**Filtek Supreme Ultra Flowable Restorative**  
 Choose \_\_\_syringe (6032) or \_\_\_capsule (6033):  
 \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4 \_\_\_B1 \_\_\_B2 \_\_\_C2 \_\_\_D2 \_\_\_W \_\_\_XW \_\_\_OA3  
**Filtek Supreme Ultra Universal Restorative**  
 Choose \_\_\_syringe (6028) or \_\_\_capsule (6029)\*\*:  
 Dentin Shades: \_\_\_A1D \_\_\_A2D \_\_\_A3D \_\_\_A4D \_\_\_B3D \_\_\_C4D \_\_\_WD  
 Body Shades: \_\_\_A1B \_\_\_A2B\*\* \_\_\_A3B \_\_\_A3.5B \_\_\_A4B \_\_\_A6B \_\_\_B1B \_\_\_B2B \_\_\_B3B \_\_\_B5B \_\_\_C1B \_\_\_C2B \_\_\_C3B \_\_\_D2B \_\_\_D3B \_\_\_WB \_\_\_XWB  
 Enamel Shades: \_\_\_A1E \_\_\_A2E \_\_\_A3E \_\_\_B1E \_\_\_B2E \_\_\_D2E \_\_\_WE \_\_\_XWE  
 Translucent Shades: \_\_\_CT \_\_\_BT \_\_\_GT \_\_\_AT  
**Filtek Universal Restorative**  
 Choose \_\_\_syringe (6555) or \_\_\_capsule (6550):  
 \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4 \_\_\_B1 \_\_\_B2 \_\_\_D3 \_\_\_XW \_\_\_PO

**F. Get one Paradigm VPS Impression Material refill**  
 \_\_\_Heavy Body (5313) \_\_\_Heavy Body Fast Set (5513)  
 \_\_\_Light Body (5311) \_\_\_Light Body Fast Set (5511)  
 \_\_\_Monophase (5312) \_\_\_Regular Body (5315)  
 \_\_\_Regular Body Fast (5515)

**G. Get one Penta Impression Material\***  
**Impregum Penta Super Quick Refill:** \_\_\_Super Quick HB (69406) \_\_\_Super Quick MB (69385)  
**Impregum Penta Soft Refill:** \_\_\_Quick Step MB (31765) \_\_\_Quick Step HB (31769) \_\_\_LB (31746) \_\_\_MB (31731) \_\_\_HB (31741)  
**Intro Kit:** \_\_\_HB/LB (P31739) \_\_\_Quick Step HB/LB (P31779) \_\_\_MB (P31735) \_\_\_Quick Step MB (P31771)  
**Impregum Penta Refill:** \_\_\_MB (31642) **Intro Kit:** \_\_\_MB (P31685)  
**Imprint 4 Penta Refill:** \_\_\_Putty (71486) \_\_\_HB (71484)\*\* \_\_\_Super Quick HB (71485) \_\_\_HB (71494)\*\*  
**Intro Kit:** \_\_\_Putty (71536) \_\_\_HB (71494)\*\* \_\_\_Super Quick HB (71495)  
**Imprint 4 Preliminary Penta Refill:** \_\_\_Regular (71521) \_\_\_Super Quick (71522)  
 Value Pack: \_\_\_Regular (71523) \_\_\_Super Quick (71524)

**H. Get one Garant Cartridge Impression Material\***  
**Impregum Super Quick Refill:** \_\_\_Super Quick LB (69379)  
**Impregum Soft Refill:** \_\_\_Tray (31785) \_\_\_Quick Step Tray (31788) \_\_\_LB (31752) \_\_\_Quick Step LB (31767)  
**Intro Kit:** \_\_\_Tray (31786) \_\_\_Quick Step Tray (31789)  
**Imprint 4 Refill:** \_\_\_HB (71492)\*\* \_\_\_LB (71488) \_\_\_RB (71487) \_\_\_Super Quick HB (71493) \_\_\_Super Quick LB (71490) \_\_\_Super Quick RB (71489)  
**Intro Kit:** \_\_\_HB (71534) \_\_\_Super Quick HB (71535)  
**Economy Pack:** \_\_\_HB (71531) \_\_\_Super Quick HB (71532)  
**Imprint 4 Preliminary Refill:** \_\_\_Regular (71539) \_\_\_Super Quick (71540)  
**Bulk Pack:** \_\_\_Regular (71548) \_\_\_Super Quick (71549)

**I. \_\_\_Get one Imprint 4 Bite VPS Bite Registration Material refill (71529)**

**J. Get one Case of Clinpro 5000 Toothpaste**  
 \_\_\_Vanilla Mint (12115)\*\* \_\_\_Spearmint (12115SM)  
 \_\_\_Bubble Gum (12115BG)

**K. Get one of the Same Varnish count-packs purchased Vanish Varnish**  
 50-Pack: \_\_\_Cherry (12149C) \_\_\_Melon (12149L)\*\* \_\_\_Mint (12149M)  
 100-Pack: \_\_\_Cherry (12150C) \_\_\_Melon (12150L)\*\* \_\_\_Mint (12150M)  
 1,000-Pack: \_\_\_Cherry (12151C) \_\_\_Melon (12151L)\*\* \_\_\_Mint (12151M) \_\_\_Cherry/Melon (12151CL) \_\_\_Cherry/Melon/Mint (12151X)  
**Fast Release Varnish**  
 50-Pack: \_\_\_Mint (12157M) \_\_\_Strawberry (12157S) \_\_\_Caramel (12157C)  
 200-Pack: \_\_\_Mint (12158M) \_\_\_Strawberry (12158S) \_\_\_Caramel (12158C)

**L. \_\_\_Get two boxes Xerostomia Relief Spray (12137)**

\*Bonus good must be of equal or lesser value of lowest price of purchased products.  
 \*\*If bonus goods associated with purchased offer(s) are not selected or noted clearly, the products indicated with asterisks (\*\*) will be automatically sent as bonus goods.

Please check the appropriate box if you do not wish to receive communications from 3M Oral Care via  fax or  email. This information will be used only if you provide information about 3M products and in accordance with 3M's privacy policy found at www.3M.com/privacy/jhtml. Offers are valid in the United States. Please allow 3–5 weeks for delivery. 3M Oral Care reserves the right to discontinue or change these offers at any time. It is the doctor's obligation to properly report and appropriately reflect any bonus product on their submissions to Medicare/Medicaid/private insurance.

3M, ESPE, Aplicap, Clicker, Clinpro, Filtek, Garant, Impregum, Imprint, Maxicap, Paradigm, Penta, RelyX, Scotchbond and Vanish are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2020. All rights reserved.

