



## Victory Series™ Active Self-Ligating Brackets

**APC™ II**  
Adhesive Coated Appliance System



### Parts List – 0.022 slot Roth\* Prescription

Maxillary			IN/OUT		0.022 in.	
Tooth	Torque	Angulation	in.	mm	L	R
Central	+12	+5	0.040	1.02	025-301	025-302
Lateral	+8	+9	0.049	1.24	025-303	025-304
Cuspid HK	-2	+10	0.032	0.81	025-307	025-308
Bicuspid HK	-7	0	0.037	0.94	025-310	025-311
Mandibular			IN/OUT		0.022 in.	
Tooth	Torque	Angulation	in.	mm	L	R
Anterior	-1	0	0.058	1.47	025-312	
Cuspid HK	-11	+7	0.030	0.76	025-315	025-316
1st Bicuspid HK	-17	0	0.039	0.99	025-319	025-320
2nd Bicuspid HK	-22	0	0.033	0.84	025-323	025-324

Adding Prefix "3" is APC™ II Adhesive System – Example: 3025-301

### Single Patient Kits



	0.022 in.
U/L 5x5 Victory Series™ Active Self-Ligating Brackets HKs	025-127
Adding Prefix "3" is APC™ II Adhesive System – Example: 3025-127	

\*3M version of this prescription. No endorsement by the Doctor is implied.