



Science.
Applied to Life.™



Reduce operating costs. Standardize on one electrode, cable and leadwire provider.

With supply chain costs at 40%¹ of a hospital's operating budget, it can be a challenge to reduce costs without compromising patient care.

Get the high quality products you expect with the 3M standardization plan:

- reduce overall usage
- reduce the number of SKUs across departments
- reduce inventory from 10-15 stocked electrodes down to 5

1. Federal Register Col 72, No. 73 Rules and Regulations (73 FR 20373)

Bring inventory down. Bring efficiency up.

No matter what department you are in, your facility can standardize on one ECG accessory provider without jeopardizing product performance. 3M has electrodes that meet many price and performance needs. And with one leadwire to trunk cable connection, facility wide, we can help make it easier and quicker for your staff to transport patients between departments.

Universal

Great for facility-wide ECG monitoring.



3M™ Foam Monitoring Electrode 2244 Radiolucent with Sticky Gel (shown above)

3M™ Foam Monitoring Electrode 2228 with Sticky Gel



3M™ Red Dot™ Monitoring Electrode 2570 Radiolucent with Abrader, Foam Tape and Sticky Gel (shown above)

3M™ Red Dot™ Monitoring Electrode 2560 with Foam Tape and Sticky Gel



3M™ Red Dot™ Repositionable Monitoring Electrode 2670 – Radiolucent (shown above)

3M™ Red Dot™ Repositionable Monitoring Electrode 2660 – Radiolucent

Resting

For EKG diagnostics typically done in the ER or EKG lab.



3M™ Red Dot™ Resting EKG Electrode 2360 Radiolucent

Pediatric

Smaller sized with a gentle adhesive.



3M™ Red Dot™ ECG Monitoring Electrode 2268-3 Radiolucent


Neonatal

Sized for newborns and infants with a gentle adhesive.

3M™ Red Dot™ Neonatal, Pre-wired, Radiolucent Soft Cloth Monitoring Electrode 2282



3M™ Red Dot™ Neonatal, Pre-Wired, Radiolucent Monitoring Electrode with Clear Tape 2269T



Disposable/Reusable

Standardize on one universal connection.

AMC&E 5-Lead Universal DIN to Pinch Disposable Leadwire with Shielded Ribbon Cable and Combiner YMDLW5P (shown below)



3M Health Care
Infection Prevention Division
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
Phone 1-800-228-3957
Web 3M.com/patientmonitoring

Ready to make a change

and start standardizing with one of each of these essential core patient monitoring categories?

Call 1-800-228-3957
or visit go.3M.com/3MRedDot

View Full Line Brochure:
go.3M.com/PatientMonitoringBrochure
View Cables and Leadwires Brochure:
go.3M.com/CablesandLeadwiresBrochure

3M and Red Dot are trademarks of 3M.
Please recycle. Printed in USA.
© 3M 2016. All rights reserved.
70-2011-5518-4