

The CRS+ Advantage: 3M™ Coding and Reimbursement System Plus

3M Coding and Reimbursement System Plus (CRS+) is a sophisticated yet easy-to-use coding solution that supports coders with powerful coding logic backed by 3M experts and nosologists. The system guides the coder through the coding process, covering the full range of situations from complex oncology treatments to repetitive or typical patient encounters.

With 3M CRS+, coders have a flexible tool that supports greater productivity—enabling direct coding for routine procedures and expert logic to help address unique or complex coding issues.

Key 3M CRS+ benefits

What it means for you



Create custom edits for greater control over the coding process at the time of coding

- Create custom edits that check coding compliance at the point-of-coding, based on your organization's unique coding needs
- Tailor edits to meet the requirements of payer contracts that are specific to your facility
- Prevent coders from completing a record until the edit is resolved (based upon an edit's severity ranking set by the organization), thus increasing accuracy and reducing rework



Use reporting tools to gain visibility into your organization's coding operations

- Increase visibility into coder performance with both inpatient and outpatient coder productivity reports
- Provide reports for ICD procedures, HCPCS procedures, DRG summaries and DX summaries
- Access reports easily without IT department intervention or support



Receive automated updates to inpatient facility Medicare reimbursement values

- Automatically receive accurate updates to reimbursement values based on Medicare regional updates
- Eliminate the need to manually update and configure your reimbursement values



Keep physician lists accurate and simple to update

- Automatically import a list of physician names from an existing system
- Save time by pre-populating physician names
- Reduce chance of identifying incorrect physician



Leverage quality indicators to ensure accurate quality metrics

- Surface AHRQ Quality Indicators and Elixhauser comorbidity values in 3M CRS+
- Integrate indicators into custom edits to ensure an organization's quality needs are managed

Call today

For more information on how 3M software and services can assist your organization, contact your 3M sales representative, call us any time toll-free at **800-367-2447**, or visit **www.3M.com/his**.